

Eye tracking capabilities in evaluating human factors during cognitive and motor tasks

INTRODUCTION

In order to test any process, product, or system, it is important to understand the underlying human factors at play. A reliable way to do this is by including biosensors in research. Simple sensors like eye tracking, facial expression analysis, galvanic skin response, and heart rate[1-4] can provide valuable information about how people work with information, for example by analyzing how they pay attention to instructions, analyze data, and identify errors. This study compares and contrasts one biosensor - eye tracking[5] - in the value it can provide for analysis of different systems.

APPARATUS & ALGORITHMS

Smart Eye AI-X (60Hz) was used for the **SBET** condition, Pupil Invisible Glasses (200 Hz) for the **Glasses** condition, and iMotions Webcam based eye tracking algorithm[6] for the **WebET** condition.




iMotions 9.3.10 was used for fixation classification, gaze mapping and AOI analysis. Fixation classification was conducted using the iMotions I-VT filter for the **SBET** and **Glasses** condition and the iMotions I-HMM filter for the **WebET** Condition. Affectiva's AFFDEX algorithm[7] was used for facial emotion detection.

EASE OF USE & DATA ACCURACY

SBET vs **WebET** showed an accuracy - scalability tradeoff. While **SBET** required a resource intensive lab set-up, it also provided accurate data. The **WebET** provided ease of scalability, but much lowered accuracy.

SBET vs **Glasses** showed an ecological validity to study resources tradeoff. The **Glasses** allowed participants to use physical blocks lending higher ecological validity, but the analysis of the individualized, unique datasets proved labor intensive.


WebET vs **Glasses** showed both to have higher ecological validity. While the **Glasses** allow participants to use actual blocks, **WebET** allows participants to conduct screen based studies from the comfort of their homes. The **Glasses** condition is resource intensive, but the **WebET** data is less accurate.



COMPARABILITY OF DATA

Area of Interest (AOI) analysis showed comparable aggregate Dwell time (~47% of time) on the building board for the **Glasses** and **SBET** condition. The **WebET** condition had a much smaller Dwell time (~22%). The heatmap also reflects the variability of the data inherent to **WebET**, suggesting that many of the true fixations on the building board may have been too noisy and outside the scope of the AOI.

FACIAL EXPRESSION ANALYSIS



Metrics, such as smile and brow furrow, measured well in all three setups. However, people had the highest percentage of facial engagement in the **Glasses** condition, where the experimenter was interacting with the participants.

People were also more likely to smile when the experimenter was present, and more likely to furrow their brow (likely a sign of concentration) when doing the task online, in a lab, or at home in the **SBET** and **WebET** conditions.

SUMMARY OF CONDITIONS COMPARING SETUPS

	SBET	Glasses	WebET
Accuracy	High accuracy	Lower accuracy, larger field of view	Lowest accuracy
Granularity of analysis	Smaller AOIs (e.g. on blocks) are possible	Increased analysis time to track all AOIs dynamically	Recommended with larger AOIs to account for accuracy errors
Ecological validity	Controlled environment	High ecological validity	Higher ecological validity, restricted to a screen
Scalability	Time consuming to collect data	Time consuming to collect and analyse	Easily scalable
Facial Expression	Integrates easily	Affected by movements	Integrates easily
Multi-sensor integration	Possible to integrate with GSR, EMG, ECG, and EEG	Easy to integrate with GSR, ECG, EMG Movement artifacts for EEG	Not possible with the Online platform
Cost effectiveness	Requires dedicated hardware and study specific resources	Requires dedicated hardware and study specific resources	Webcam based. Minimum human resources.

READ THE FULL ARTICLE HERE



REFERENCES

1. Mele & Federici, *Cogn Process* (2012)
2. Hazlett, *CHI '06* (2006)
3. Hardy et al., *Proc Hum Factors Ergo Soc*, (2013)
4. Hankins & Wilson, *Aviat Space Environ Med* (1998)
5. Holmqvist et al., *Behavior research method* (2022)
6. iMotions, "Webcam ET Accuracy Whitepaper v2" (2022)
7. Bishay et al., *arXiv.org* (2022)