




Human Factors

in Automotive Human-Machine Interface
(HMI) Design



Powering Human Insights

Confidential © 2024 Redistribution is not permitted without written permission from IMotions

Introduction

Human-Machine Interface (HMI) design is a crucial aspect of the modern automotive design process. Through the use of biosensors, iMotions' emotion analytics and Smart Eye's eye-tracking technology redefine driver experience optimization. Learn how these advancements ensure interfaces are intuitive, safe, and tailored to human behavior for a seamless vehicle-user interaction.

In automotive design, the integration of Human-Machine Interface (HMI) systems stands as a rapidly evolving testament to the evolution of vehicle-user interaction. The complexity of these systems necessitates a nuanced understanding of human factors engineering, ensuring that the driver's experience is intuitive, safe, and efficient.

Companies like iMotions and its sister company Smart Eye are at the forefront of this innovation, leveraging their expertise in behavioral research and eye-tracking technologies to redefine HMI design principles.


iMotions Quick Facts

iMotions is the leading software company for multimodal biometric research solutions.

Using biometric data in automotive research allows for gathering quantitative data together with qualitative data.

- Monitor and record physiological state
- Assess cognitive load during different tasks
- Multiple camera feed synchronization with annotation
- Visualization and analysis of eye tracking data
- API for realtime data loops

Collect, synchronize and analyze multiple data sources with iMotions



The Role of Human Factors in HMI Design

Human factors in HMI design refer to the consideration of human capabilities and limitations in the context of interacting with automotive user interfaces, which essentially is everything the driver interacts with while driving.

This involves a delicate balance between ergonomic design, cognitive psychology, and user experience (UX) research to create interfaces that are not only user-friendly but also augment safety and performance all the while making the car intuitive and pleasant to operate. The goal is to design systems that align with natural human behaviors and cognitive processes to minimize the potential for errors and maximize usability. But what are the different aspects of Human factors when applied to an HMI design process?

Cognitive Load and Information Processing

One of the critical considerations in automotive HMI design is managing cognitive load in the driver — the amount of mental effort being used in the working memory. Human factors experts strive to design interfaces that facilitate easy information processing, minimizing cognitive load to prevent overload and confusion in order to prevent dangerous situations from arising.

Modern steering wheels doubles as control panels to non-essential operations for the driver so as to increase safety and comfort.



This involves organizing information logically, avoiding unnecessary complexity, and leveraging visual and interactive hierarchies to prioritize critical functions. A good example of this is the various buttons on a steering wheel controlling a number of non-driving-critical operations, such as AC, music, cruise control, so that drivers can keep their eyes on the road while engaging with these features.

Ergonomics and Physical Interaction

Ergonomics plays a vital role in automotive HMI design, focusing on the physical interaction between the user and the interior of a car. This includes the placement of controls and displays, the physical effort required to operate them, and the feedback they provide. Good ergonomic design ensures that the interface can be used comfortably and efficiently, reducing the risk of physical strain, error or accidents.

User-Centered Design and Usability

At the heart of human factors is the principle of user-centered design — the process of involving users throughout the design and development phases to create more effective interfaces. This approach emphasizes usability, aiming to make systems easy to learn, efficient to use, and satisfying. User-centered design involves iterative testing with real users, allowing designers to refine the interface based on direct feedback and observations of user behavior.

Emotional and Psychological Considerations

Beyond cognitive and physical aspects, human factors in HMI design also consider the emotional and psychological impact of interfaces on users. This includes understanding how design elements can affect mood and stress levels, as well as designing interfaces that foster positive emotional connections with the user. For example, the use of color, shape, and texture can influence a user's perception and emotional response to a system.

Safety and Accessibility

Safety is a paramount concern in HMI design, with human factors expertise being crucial in minimizing errors and accidents. This involves designing interfaces that are clear, predictable, and forgiving of user mistakes. Additionally, accessibility is a key aspect of human factors, ensuring that systems can be used by people with a wide range of abilities, including those with disabilities.

Emotion Analytics in HMI Design – The iMotions Way

As you might know, iMotions specializes in integrating emotion analytics into the automotive HMI design process. Their approach centers on understanding the emotional and cognitive states of drivers in real-time, leveraging advanced sensor technologies to capture data on eye movement, facial expressions, and physiological responses. This data-driven approach allows designers to identify how drivers interact with different HMI elements, pinpointing areas of cognitive overload, distraction, or stress.


By utilizing iMotions' technology, automotive designers can make informed decisions about the layout, complexity, and functionality of HMI systems. For example, insights into gaze patterns can inform the optimal placement of critical information on dashboards, ensuring that drivers can access the information they need without

diverting attention from the road. Similarly, monitoring physiological responses during interactions with infotainment systems can help in designing interfaces that minimize cognitive strain.

iMotions in Automotive Simulators

Integrating iMotions with automotive simulators marks a revolutionary step in enhancing the Human-Machine Interface (HMI) design process. This synergy allows for the creation of immersive, controlled environments where designers can experiment with and evaluate interface layouts, functionalities, and user interactions under realistic conditions.

By leveraging iMotions' advanced emotion analytics and biometric sensing



technologies within automotive simulators, researchers gain another layer of insights into the driver's cognitive load, emotional states, and physiological responses to various HMI designs.

The application of iMotions Software in this context can involve a diverse range of biometric measurements such as; tracking eye movements, facial expressions, heart rate variability, and EEG readings among other metrics, to ascertain how drivers interact with and react to different interface elements during simulated driving tasks.


This real-time data provides a rich, nuanced understanding of driver behavior, enabling designers to identify which aspects of the HMI design facilitate intuitive use and which may cause confusion or distraction.

Furthermore, the feedback loop created by integrating iMotions with automotive simulators accelerates the HMI design iteration process. Designers can quickly prototype and test modifications, effectively using empirical evidence to guide the development of interfaces that are not only aesthetically pleasing but also optimize safety and user satisfaction.

This approach not only enhances the design process but also contributes to the creation of automotive interfaces that are more aligned with human capabilities and limitations, paving the way for safer, more intuitive driving experiences.

Why Biometrics in Automotive

Biometrics allow researchers to take a multi-dimensional approach to understanding and explaining human behavior.



Smart Eye Pro: Pioneering Eye-Tracking Technology in Automotive HMI Design

Smart Eye's contribution to HMI design lies in its advanced eye-tracking technology through their Smart Eye Pro system. The multi-camera system allows HMI designers and researchers to precisely track where and how long a driver looks at different parts of the vehicle interior, and how they engage with the cabin.

Smart Eye Pro provides invaluable insights into attention distribution and visual engagement. This technology is instrumental in developing adaptive HMI systems that can dynamically adjust based on the driver's focus, enhancing both safety and user experience.



Example of setup with iMotions and Smart Eye Pro

Synergizing Strengths for Future HMI Innovations


The strengths of iMotions and Smart Eye in emotion analytics and eye-tracking technology, respectively, offer the potential for synergy for the future of automotive HMI design. By combining these technologies, designers can gain a comprehensive understanding of the driver's physical and emotional state, leading to the creation of HMIs that are not only more intuitive and engaging but also contribute to safer driving environments.

The Path Forward in HCI Design

As we move towards a future where vehicles become increasingly autonomous, the role of HMI design in ensuring seamless interaction between humans and machines becomes ever more critical. Technologies such as the ones provided

by iMotions and Smart Eye, not only help in optimizing current designs but also lay the groundwork for the next generation of automotive interfaces, promising a future where vehicles are not just modes of transport but partners in enhancing human capability and safety on the road.

Easy Video Observation with iMotions Lab



Synchronize and annotate multiple video sources inside iMotions Lab, allowing for behavioral coding and general observations to be recorded with ease.




The quick brown ...

Add additional data layers such as automatic transcription of spoken words with sentiment analysis for even greater insights and accurate recording of observations.

iMotions Software Solution

Multimodal research in any environment

iMotions reduces the complexity of carrying out multimodal research, enabling a wide array of sensors to be seamlessly connected. By combining these different physiological measurements, it's possible to get a better understanding of human behavior in any environment.



iMotions enables multimodal research to be carried out in an array of research scenarios.

Lab-based studies




Natural environments



VR environments



iMotions supports leading 3rd party sensor products. Additional sensors can be integrated via our API.



Biosensors

iMotions has a suite of partners that provide biosensor hardware like EEG, electrodermal activity (EDA), ECG, EMG, etc that are suitable for automotive research.

- Eye tracking: SMI, Smart Eye, EyeTech, Gazepoint
- EEG: Advanced Brain Monitoring (ABM), NeuroElectrics, Brain Products, Emotiv, OpenBCI
- EDA: BIOPAC, Shimmer, Empatica
- ECG: BIOPAC, Shimmer
- EMG: BIOPAC, Shimmer
- Respiration: BIOPAC

Additional sensor integration

iMotions has a powerful API equipping researchers with the tools to:

- Integrate new sensors with real-time data capture in iMotions
- Live synchronize all data streams
- Live forward all synchronized data streams
- Control the iMotions application via remote control
- Supports standard sensor protocols like LSL and TTL

EEG



EDA



ECG / EMG



David Hurwitz, Director of Kiewit Center for Infrastructure and Transportation Research at Oregon State University, describes his experience setting up a cutting edge multimodal transportation simulation lab with iMotions.



Scan to watch



Oregon State
University

Want to know more?

GET IN TOUCH

IMOTIONS®

Copenhagen, Denmark

Kristen Bernikows Gade 6
4th floor
København K, 1105
TEL +45 71 998 098

Boston, USA

38 Chauncy Street
Floor 8, Suite 800
Boston, MA 02111
TEL +1 617-520-4958

Synchronize, Visualize and Analyze your research in Eye Tracking, Facial Expression Analysis, Galvanic Skin Response, Surveys, EEG and much more in one software platform.

www.imotions.com

China

NO.1 Fortune Avenue,
Room 2902
Yubei District, Chongqing
TEL +886 931684806

